

Screen, Select, Enhance: An Integrated Workflow for ADCC-Enhanced Antibody Development and Validation

evitria
your proteins, reliably delivered


Julia Pizzolato¹, Julian Wichmann², Tobias Luttermann² | ¹ evitria AG, Brunastrasse 185, CH-8951 Fahrweid, Switzerland
² Biofidus AG, Guetersloher Strasse 54, DE-33649 Bielefeld, Germany

From high-throughput screening to validated ADCC-enhanced lead candidates

1 Introduction

This study focuses on evitria's high-throughput (HTP) antibody expression pipeline from candidate selection and scale-up of ADCC-enhanced variants to mechanistic characterization by our partner Biofidus. Leveraging robust, standardized processes with high flexibility, evitria offers expertise in bispecific development, high-throughput expression, and Fc-engineered antibodies. We also in-license technologies such as ProBioGen's GlymaxX™ system, which afucosylates the antibody Fc to enhance ADCC effector function.

Biofidus AG, a leading (bio)analytical laboratory, provides comprehensive bioanalytical services as a dynamic one-stop partner for biologics. As a member of the Midas Pharma Group, we are embedded in a global network, expanding our capabilities and international reach to deliver even greater value to our partners. At Biofidus, we use our Biacore 1K surface plasmon resonance (SPR) platform to deliver fast, label-free binding data for an in-depth characterization of biologics.




2 GlymaxX™ technology


GlymaxX™ redirects fucose synthesis using a heterologous enzyme that depletes the intracellular fucose pool. This enables production of both fucosylated and afucosylated antibodies from the same cell line, with no amino acid sequence change.

Afucosylated antibodies show up to 40-fold increased affinity to FcγRIIIA/B, resulting in up to 100x higher ADCC activity. A unique glycan at Asn162 in FcγRIIIA/B (not present in other FcγR subtypes) collides with the core fucose on IgG-Fc.

Afucosylation therefore specifically boosts ADCC via FcγRIIIA-expressing NK cells, without altering other receptor interactions.




3 End-to-end workflow from HTP protein expression and screening to functional data with evitria and Biofidus




4 Selection of top candidates and expression +/- GlymaxX™

After initial HTP expression of a 48-constructs panel of different binder variants (Figure A), the best expressing and most stable construct was selected and up-scaled to 250 mL expression (+/- GlymaxX™, +/- LALAPG) followed by protein A purification (Table B)




→ High expression levels achieved across constructs without need for polishing
→ Final candidate rapidly identified for downstream characterization


5 Validating target binding of final lead candidates via SPR




6 GlymaxX™ leads to increased binding to FcγRIIIA-176V



7 GlymaxX™ does not affect binding to FcγRI



8 GlymaxX™ effect can be recapitulated with mouse IgG2a



9 Conclusions

- ✓ evitria's integrated HTP-to-lead pipeline compresses ADCC antibody development timelines without sacrificing analytical rigor by screening constructs directly in CHO
- ✓ GlymaxX™-mediated afucosylation selectively increases FcγRIIIA binding (4.4-fold K_D decrease) without affecting FcγRI binding or target engagement
- ✓ LALAPG and GlymaxX are orthogonal, stackable Fc engineering tools: their opposing effects on FcγRIIIA binding allow fine-tuned effector function across a wide dynamic range
- ✓ Comparable results in a mouse background confirm cross-species applicability of GlymaxX™, de-risking translational studies
- ✓ evitria-Biofidus partnership closes the loop from sequence to functional ADCC data under a single coordinated workflow

10 Material & Methods

Protein expression and purification
cDNA was cloned into evitria's vector system. Suspension-adapted CHO K1 cells were transfected using eviFect in serum-free media. Supernatants were harvested by centrifugation, and antibodies were purified using MabSelect™ SuRe™ Prot A columns.

SPR spectroscopy
Samples were analyzed on a Biacore 1K system (Cytiva). His-tagged human FcγRIIIA-176V, target antigen, or FcγRI were captured on a CM5 chip via immobilized anti-histidine antibody. Double-referenced sensorgrams were evaluated using Biacore Insight Evaluation Software with a 1:1 binding model for kinetic constants and steady-state analysis for equilibrium dissociation constants (K_D).